**Breeze Nasal** : Referentie : Maat S= M-268020-10 M= M-268030-10 L= M-268040-10

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**RECOMMENDED USE**

The Breeze Nasal Mask is designed for nasal ventilation via the use of a barometric device delivering Continuous Positive Airway Pressure (CPAP) or bi-level pressure (BiPAP).

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
Store the mask in a clean and dry area, out of the sun. To protect the environment, the mask must be disposed of in a waste bin. Do not discard improperly. Storage life : 5 years from the date of manufacture.

**TECHNICAL CHARACTERISTICS**










Treatment pressure :	3 cmH <sub>2</sub> O to 20 cmH <sub>2</sub> O
Operating temperature :	5° C to +35° C
Storage temperature :	-20° C to +60° C
Pressure drop measurement :	at 50 l/min 0.28 cmH <sub>2</sub> O and at 100 l/min 1.09 cmH <sub>2</sub> O
Noise level :	28.0 dBA to 10cmH <sub>2</sub> O
Materials used in contact with the patient : silicone, Oeko-Tex textile and SBR elastomer	

**PUTTING THE MASK IN PLACE**

Prepare the harness by detaching the Velcro at A1 and A2. (fig. 1). Insert the A1 Velcro strip in B1 and close (fig. 3 and 4). Repeat for A2 and B2 (fig. 5 and 6). Place the mask on your head and adjust the tightness of the upper tethers of the harness (fig. 7). Clip the A3 loop to the B3 hook of the frame. Repeat for A4 and B4. (fig. 8). Adjust the tightness of the A1, A2, A3 and A4 Velcros to position the harness properly. (fig. 9 and 10) Make sure that the back of the harness is properly positioned, to prevent slipping. In case of daily use, just clip the A3 loop to the B3 hook of the frame, or A4 to B4 (fig.8), and unhook it at the end of treatment.

List of components:	
1 Silicone cup	5 Plug
2 Nasal frame	6 Adapter
3 Forehead support	7 Tether loop
4 Elbow	8 Harness

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**RECOMMENDED USE**

The Breeze Nasal + Mask is designed for nasal ventilation via the use of a barometric device delivering Continuous Positive Airway Pressure (CPAP) or bi-level pressure (BiPAP).

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**WARNINGS**

- Do not block or attempt to obstruct the leak ports.
- This mask can only be used with CPAP or BiPAP devices recommended by a physician or lung specialist.
- The mask cannot be used if the device has not been started up and is functioning normally.
- Explanation:** CPAP and BiPAP devices must be used with specific masks equipped with leak ports enabling the continuous discharge of exhaled air. When the device is functioning normally, fresh air from the device will evacuate the exhaled air through the leak ports of the mask. However, when the device is not operational, there is not enough fresh air coming into the mask, and the exhaled air could be inhaled again. Breathing in exhaled air for longer than a few minutes can, in certain circumstances, lead to suffocation. This applies to most CPAP or BiPAP devices.
- Stop using this mask if you experience any adverse reaction whatsoever during usage (claustrophobia, irritation, pain, conjunctivitis, etc...).
- Should this occur, consult your GP.
- Using a mask can cause pain in the gums, the teeth or the jaw or can aggravate dental problems. If this occurs, please see a dentist.
- Do not use the mask if the packaging has been opened or damaged. Replace the mask or the equipment if it is damaged or distorted.
- Clean the mask before using it for the first time.
- Do not use if the patient cannot remove the mask unaided.
- If oxygen is to be introduced, follow the instructions and safety precautions associated with the use of oxygen.
- To minimize leakage, different cup and forehead support sizes are available. Ask your health care provider.

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**STORAGE**


The mask must be washed and dried before being put away.

Store the mask in a clean and dry area, out of the sun.

To protect the environment, the mask must be disposed of in a waste bin. Do not discard improperly.  
 Storage life : 5 years from the date of manufacture.

**TECHNICAL CHARACTERISTICS**

Treatment pressure :	3 cmH <sub>2</sub> O to 20 cmH <sub>2</sub> O
Operating temperature :	5° C to +35° C
Storage temperature :	-20° C to +60° C
Pressure drop measurement :	at 50 l/min 0.24 cmH <sub>2</sub> O and at 100 l/min 1.04 cmH <sub>2</sub> O
Noise level :	28.0 dBA to 10cmH <sub>2</sub> O
Materials used in contact with the patient :	silicone, Oeko-TEX textile and SBR elastomer


Pressure-flow curve :

3 cmH <sub>2</sub> O :	15.4 l/Min.
7.25 cmH <sub>2</sub> O :	25.2 l/Min.
11.5 cmH <sub>2</sub> O :	33.1 l/Min.
15.75 cmH <sub>2</sub> O :	42.3 l/Min.
20 cmH <sub>2</sub> O :	52.6 l/Min.

**PUTTING THE MASK IN PLACE**

Prepare the harness by detaching the Velcro at A1 and A2. (fig. 1).  
 Insert the A1 Velcro strip in B1 and close (fig. 3 and 4). Repeat for A2 and B2 (fig. 5 and 6).  
 Place the mask on your head and adjust the tightness of the upper tethers of the harness (fig. 7).  
 Clip the A3 loop to the B3 hook of the frame. Repeat for A4 and B4 (fig. 8).  
 Adjust the tightness of the A1, A2, A3 and A4 Velcros to position the harness properly. (fig. 9 and 10)  
 Make sure that the back of the harness is properly positioned, to prevent slipping.  
 In case of daily use, just clip the A3 loop to the B3 hook of the frame, or A4 to B4 (fig. 8), and unhook it at the end of treatment.

List of components :	
1 Silicone cup	6 Adapter
2 Nasal frame	7 Tether loop
3 Forehead support size M	8 Harness
4 Elbow	9 Forehead support size S
5 Plug	10 Forehead support size L


**Breeze Nasal +**

www.sefam-medical.com

SEFAM-Customer Service  
 10, allée Pelletier Doisy  
 54600 Villers-lès-Nancy  
 France

Tél : +33 383 44 85 00  
 Fax : +33 383 44 85 01

Email :  
 customerservice@sefam-medical.com

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**RECOMMENDED USE**

The Breeze Nasal Comfort Mask is designed for nasal ventilation via the use of a barometric device delivering Continuous Positive Airway Pressure (CPAP) or bi-level pressure (BiPAP).  
The Breeze Nasal Comfort Mask must be:

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**WARNINGS**

- Do not block or attempt to obstruct the leak ports.
  - This mask can only be used with CPAP or BiPAP devices recommended by a physician or lung specialist.
  - The mask cannot be used if the device has not been started up and is functioning normally.
- Explanation: CPAP and BiPAP devices must be used with specific masks equipped with leak ports enabling the continuous discharge of exhaled air. When the device is functioning normally, fresh air from the device will evacuate the exhaled air through the leak ports of the mask. However, when the device is not operational, there is not enough fresh air coming into the mask, and the exhaled air could be inhaled again. Breathing in exhaled air for longer than a few minutes can, in certain circumstances, lead to suffocation. This applies to most CPAP or BiPAP devices.*
- Stop using this mask if you experience any adverse reaction whatsoever during usage (claustrophobia, irritation, pain, conjunctivitis, etc...).
  - Should this occur, consult your GP.
  - Using a mask can cause pain in the gums, the teeth or the jaw or can aggravate dental problems. If this occurs, please see a dentist.
  - Do not use the mask if the packaging has been opened or damaged. Replace the mask or the equipment if it is damaged or distorted.
  - Clean the mask before using it for the first time.
  - Do not use if the patient cannot remove the mask unaided.
  - If oxygen is to be introduced, follow the instructions and safety precautions associated with the use of oxygen.
  - To minimize leakage, different cup and forehead support sizes are available. Ask your health care provider.

**MAINTENANCE / CLEANING**

- The mask must be cleaned every day in accordance with the following instructions:
- Separate the mask from any other accessory (harness, corrugated tube).
  - Wash the different parts of the mask with warm water and a mild detergent (5 <pH < 8) such as soap or washing-up liquid.
  - Rinse every part with clean water.
  - Allow the equipment to dry completely on a clean cloth out of the sun before using it again.
  - Wash the harness by hand once a week.


- Note:**
- Do not use bleach, alcohol or cleaning solutions containing alcohol, shampoo or strongly scented products.
  - If part of the system becomes visibly deteriorated (crack, fissure, tear, etc...), the component must be replaced.

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**TECHNICAL CHARACTERISTICS**












Treatment pressure :	3 cmH <sub>2</sub> O to 20 cmH <sub>2</sub> O
Operating temperature :	5° C to +35° C
Storage temperature :	-20° C to +60° C
Pressure drop measurement :	at 50 l/min 0.2 cmH <sub>2</sub> O and at 100 l/min 0.85 cmH <sub>2</sub> O
Noise level :	27.5 dBA to 10cmH <sub>2</sub> O
Materials used in contact with the patient : silicone, Oeko-Tex textile and SBR elastomer	

**PUTTING THE MASK IN PLACE**

Prepare the harness by detaching the Velcro at A1 and A2.  
Insert the A1 Velcro strip in B1 and close (fig. 3 and 4). Repeat for A2 and B2 (fig. 5 and 6).  
Place the mask on your head and adjust the tightness of the upper tethers of the harness (fig. 7).  
Clip the A3 tether to the B3 hook of the frame. Repeat for A4 and B4. (fig. 8).  
Adjust the tightness of the A1, A2, A3 and A4 Velcros to position the harness properly. (fig. 9 and 10)  
Make sure that the back of the harness is properly positioned, to prevent slipping.  
In case of daily use, just clip the A3 tether to the B3 hook of the frame, or A4 to B4 (fig.8), and remove it at the end of treatment.

List of components :	
1 Silicone cup	6 Plug
2 Internal frame	7 Tether clip
3 External frame	8 Adapter
4 Forehead support	9 Harness
5 Elbow	

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
La mascarilla Breeze Nasal Comfort debe ser:

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Please read these instructions carefully before using your mask.

**RECOMMENDED USE**

The Breeze Facial + Mask is designed for non invasive ventilation via the use of a barometric device delivering Continuous Positive Airway Pressure (CPAP) or bi-level pressure (BiPAP).

The Breeze Facial + Mask must be:

- used by adult patients (>30Kg) for whom a positive airway pressure treatment was prescribed.
- used by a single patient treated at home or in a health care facility.

**WARNINGS**

- Do not block or attempt to obstruct the leak ports.
  - This mask can only be used with CPAP or BiPAP devices recommended by a physician or lung specialist.
  - The mask cannot be used if the device has not been started up and is functioning normally.
- Explanation: CPAP and BiPAP devices must be used with specific masks equipped with leak ports enabling the continuous discharge of exhaled air. When the device is functioning normally, fresh air from the device will evacuate the exhaled air through the leak ports of the mask. However, when the device is not operational, there is not enough fresh air coming into the mask, and the exhaled air could be inhaled again. Breathing in exhaled air for longer than a few minutes can, in certain circumstances, lead to suffocation. This applies to most CPAP or BiPAP devices.*
- Breeze Facial Mask + also has a safety expiratory valve (AAV Valve) which allows you to inspire ambient air in case of failure of the treatment device. Do not block or remove the valve.
  - Stop using this mask if you experience any adverse reaction whatsoever during usage (claustrophobia, irritation, pain, conjunctivitis, etc...). Should this occur, consult your GP.
  - Using a mask can cause pain in the gums, the teeth or the jaw or can aggravate dental problems. If this occurs, please see a dentist.
  - Do not use the mask if the packaging has been opened or damaged. Replace the mask or the equipment if it is damaged or distorted.
  - Clean the mask before using it for the first time.
  - Do not use if the patient cannot remove the mask unaided.
  - If oxygen is to be introduced, follow the instructions and safety precautions associated with the use of oxygen.
  - To minimize leakage, different cup and forehead support sizes are available. Ask your health care provider.
  - Do not use mask in case of nausea or vomiting.

**MAINTENANCE / CLEANING**

The mask must be cleaned every day in accordance with the following instructions:

- Separate the mask from any other accessory (harness, corrugated tube).
- Wash the different parts of the mask with warm water and a mild detergent (5 <pH < 8) such as soap or washing-up liquid.
- Rinse every part with clean water.
- Allow the equipment to dry completely on a clean cloth out of the sun before using it again.
- Wash the harness by hand once a week.

**Note:**


- Do not use bleach, alcohol or cleaning solutions containing alcohol, shampoo or strongly scented products.
- If part of the system becomes visibly deteriorated (crack, fissure, tear, etc...), the component must be replaced.

**STORAGE**

The mask must be washed and dried before being put away. Store the mask in a clean and dry area, out of the sun. To protect the environment, the mask must be disposed of in a waste bin. Do not discard improperly. Storage life : 5 years from the date of manufacture.

**TECHNICAL CHARACTERISTICS**

Treatment pressure :	3 cmH <sub>2</sub> O to 20 cmH <sub>2</sub> O
Operating temperature :	5° C to +35° C
Storage temperature :	-20° C to +60° C
Pressure drop measurement :	at 50 l/min 0,28 cmH <sub>2</sub> O and at 100 l/min 1,09 cmH <sub>2</sub> O
Noise level :	27,5 dBA to 10cmH <sub>2</sub> O
Materials used in contact with the patient : silicone, Oeko-Tex textile and SBR elastomer	

Pressure-flow curve :










3 cmH <sub>2</sub> O :	29.6 l/Min.
7.25 cmH <sub>2</sub> O :	47.0 l/Min.
11.5 cmH <sub>2</sub> O :	59.4 l/Min.
15.75 cmH <sub>2</sub> O :	70.5 l/Min.
20 cmH <sub>2</sub> O :	79.9 l/Min.

**PUTTING THE MASK IN PLACE**

Prepare the harness by detaching the Velcro at A1 and A2. Insert the A1 Velcro strip in B1 and close (fig. 3 and 4). Repeat for A2 and B2 (fig. 5 and 6). Place the mask on your head and adjust the tightness of the upper tethers of the harness (fig. 7). Clip the A3 tether to the B3 hook of the frame. Repeat for A4 and B4. (fig. 8). Adjust the tightness of the A1, A2, A3 and A4 Velcros to position the harness properly. (fig. 9 and 10) Make sure that the back of the harness is properly positioned, to prevent slipping. In case of daily use, just clip the A3 tether to the B3 hook of the frame, or A4 to B4 (fig.8), and remove it at the end of treatment.

List of components :

1 Silicone cup	7 AAV Valve
2 Nasal frame	8 Iether loop
3 Forehead support size M	9 Harness
4 Elbow	10 Forehead support size S
5 Plug	11 Forehead support size L
6 Adapter	

**Breeze Facial +**

www.sefam-medical.com

Tél : +33 383 44 85 00  
Fax : +33 383 44 85 01

Email :  
customerservice@sefam-medical.com

SEFAM-Customer Service  
10, allée Pelletier Doisy  
54600 Villers-lès-Nancy  
France

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**RECOMMENDED USE**  
The Breeze Facial Comfort Mask is designed for non invasive ventilation via the use of a barometric device delivering Continuous Positive Airway Pressure (CPAP) or bi-level pressure (BiPAP).  
The Breeze Facial Comfort Mask must be:

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**WARNINGS**

- Do not block or attempt to obstruct the leak ports.
- This mask can only be used with CPAP or BiPAP devices recommended by a physician or lung specialist.
- The mask cannot be used if the device has not been started up and is functioning normally.  
**Explanation:** CPAP and BiPAP devices must be used with specific masks equipped with leak ports enabling the continuous discharge of exhaled air. When the device is functioning normally, fresh air from the device will evacuate the exhaled air through the leak ports of the mask. However, when the device is not operational, there is not enough fresh air coming into the mask, and the exhaled air could be inhaled again. Breathing in exhaled air for longer than a few minutes can, in certain circumstances, lead to suffocation. This applies to most CPAP or BiPAP devices.
- Breeze Facial Mask Comfort also has a safety expiratory valve (AAV Valve) which allows you to inspire ambient air in case of failure of the treatment device. Do not block or remove the valve.
- Stop using this mask if you experience any adverse reaction whatsoever during usage (claustrophobia, irritation, pain, conjunctivitis, etc...). Should this occur, consult your GP.
- Using a mask can cause pain in the gums, the teeth or the jaw or can aggravate dental problems. If this occurs, please see a dentist.
- Do not use the mask if the packaging has been opened or damaged. Replace the mask or the equipment if it is damaged or distorted.
- Clean the mask before using it for the first time.
- Do not use it if the patient cannot remove the mask unaided.
- If oxygen is to be introduced, follow the instructions and safety precautions associated with the use of oxygen.
- To minimize leakage, different cup and forehead support sizes are available. Ask your health care provider.
- Do not use mask in case of nausea or vomiting.

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
**STORAGE**

The mask must be washed and dried before being put away.  
Store the mask in a clean and dry area, out of the sun.  
To protect the environment, the mask must be disposed of in a waste bin. Do not discard improperly.  
Storage life : 5 years from the date of manufacture.

**TECHNICAL CHARACTERISTICS**

Treatment pressure :	3 cmH <sub>2</sub> O to 20 cmH <sub>2</sub> O
Operating temperature :	5° C to +35° C
Storage temperature :	-20° C to +60° C
Pressure drop measurement :	at 50 l/min 0.28 cmH <sub>2</sub> O and at 100 l/min 1.09 cmH <sub>2</sub> O
Noise level :	25.0 dBA to 10cmH <sub>2</sub> O
Materials used in contact with the patient : silicone, Oeko-Tex textile and SBR elastomer	


**PUTTING THE MASK IN PLACE**

Prepare the harness by detaching the Velcro at A1 and A2.  
Insert the A1 Velcro strip in B1 and close (fig. 3 and 4). Repeat for A2 and B2 (fig. 5 and 6)  
Place the mask on your head and adjust the tightness of the upper tethers of the harness (fig. 7)  
Clip the A3 tether to the B3 hook of the frame. Repeat for A4 and B4. (fig. 8)  
Adjust the tightness of the A1, A2, A3 and A4 Velcros to position the harness properly. (fig. 9 and 10)  
Make sure that the back of the harness is properly positioned, to prevent slipping.  
In case of daily use, just clip the A3 tether to the B3 hook of the frame, or A4 to B4 (fig. 8), and remove it at the end of treatment.

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La mascarilla Breeze Facial Comfort debe ser:

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**RECOMMENDED USE**

The Breeze Pillows Mask is designed for nasal ventilation via the use of a barometric device delivering Continuous Positive Airway Pressure (CPAP) or bi-level pressure (BiPAP).  
The Breeze Pillows Mask must be:  
• used by adult patients (>30Kg) for whom a positive airway pressure treatment was prescribed.  
• used by a single patient treated at home or in a health care facility.

**WARNINGS**

- Do not block or attempt to obstruct the leak ports.
- This mask can only be used with CPAP or BiPAP devices recommended by a physician or lung specialist.
- The mask cannot be used if the device has not been started up and is functioning normally.

*Explanation: CPAP and BiPAP devices must be used with specific masks equipped with leak ports enabling the continuous discharge of exhaled air. When the device is functioning normally, fresh air from the device will evacuate the exhaled air through the leak ports of the mask. However, when the device is not operational, there is not enough fresh air coming into the mask, and the exhaled air could be inhaled again. Breathing in exhaled air for longer than a few minutes can, in certain circumstances, lead to suffocation. This applies to most CPAP or BiPAP devices.*

- Stop using this mask if you experience any adverse reaction whatsoever during usage (claustrophobia, irritation, pain, conjunctivitis, etc...). Should this occur, consult your GP.
- Do not use the mask if the packaging has been opened or damaged. Replace the mask or the equipment if it is damaged or distorted.
- Clean the mask before using it for the first time.
- Do not use if the patient cannot remove the mask unaided.
- If oxygen is to be introduced, follow the instructions and safety precautions associated with the use of oxygen.

**MAINTENANCE / CLEANING**

- The mask must be cleaned every day in accordance with the following instructions:
- Separate the mask from any other accessory (harness, corrugated tube).
  - Wash the different parts of the mask with warm water and a mild detergent (5 < pH < 8) such as soap or washing-up liquid.
  - Rinse every part with clean water.
  - Allow the equipment to dry completely on a clean cloth out of the sun before using it again.
  - Wash the harness by hand once a week.

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
**TECHNICAL CHARACTERISTICS**

Treatment pressure :	3 cmH <sub>2</sub> O to 40 cmH <sub>2</sub> O
Operating temperature :	5° C to +35° C
Storage temperature :	-20° C to +60° C
Pressure drop measurement :	at 50 l/min 2.1 cmH <sub>2</sub> O and at 100 l/min 8 cmH <sub>2</sub> O
Noise level :	22.5 dBA to 10cmH <sub>2</sub> O
Materials used in contact with the patient : silicone, Oeko-Tex textile and SBR elastomer	


Pressure-flow curve :	
3 cmH <sub>2</sub> O :	13.5 l/Min.
10 cmH <sub>2</sub> O :	26.5 l/Min.
20 cmH <sub>2</sub> O :	39.6 l/Min.
30 cmH <sub>2</sub> O :	50.0 l/Min.
40 cmH <sub>2</sub> O :	58.8 l/Min.

**PUTTING THE MASK IN PLACE**

Prepare the nasal pads fig. 1  
Insert the B1 Velcro strip in A1 and close (fig. 2,3 and 4). Repeat for A1\* and B1\* (fig. 5)  
Insert the C1 Velcro strip in A2 and close (fig. 6). Repeat for C1 and A2\* (fig. 7)  
Position the mask on your head and adjust the tightness at the top and rear, fig. 8, 9 and 10  
Adjust the positioning angle of the nasal pads (Fig.10) and your mask is ready for use (fig. 11).  
Install the tube over your head (fig 12). If the nasal pads pull on your nose, adjust the tube length between your nose and the upper tether (fig. 13)  
In case of daily use, elasticity of the harness allows to remove and replace the mask without readjustment.

List of components :

1 Upper harness	6 Connection ring
2 Rear harness	7 Angled connector
3 Nasal pad	8 Silicone tube
4 Nasal pad support	9 Pivot connector
5 Main frame	10 Tube connector

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
EC REP Mdi Europa GmbH  
Langenhagener Str. 71  
30855 Langenhagen  
GERMANY

HSINER Co. Ltd.  
No 312, Jongshan Rd.  
Shengang Dist.  
Taichung City 429

## PACKAGING / LABEL BREEZE FACIAL COMFORT + MASK

BREEZE FACIAL+ MASK

M-268020-20 (size S), M-268030-20 (size M), M-268040-20 (size L)

### Packaging



# PACKAGING / LABEL BREEZE FACIAL COMFORT + MASK

BREEZE FACIAL+ MASK

M-268020-20 (size S), M-268030-20 (size M), M-268040-20 (size L)

## Definition of label stuck on the plastic bag

<b>BREEZE FACIAL+ Mask</b>			
Masque facial / Facial mask / Gesichtsmaske / Gezichtsmasker / Máscara facial / Maschera facciale			
 0434			HSINER Co., Ltd. No.312, Jhongshan Rd., Shengang Dist. Taichung City 429
<b>REF</b> M-268020-20			MADE IN TAIWAN
<b>LOT</b> (10) 2013072510	<b>EC REP</b>		Mdi Europa GmbH Langenhagener Str. 71 30855 Langenhagen - GERMANY
 (11) 141017		<b>RX ONLY</b>	 -20°C / 60°C
 (17) 191017			

## PACKAGING / LABEL BREEZE FACIAL COMFORT MASK

### BREEZE FACIAL COMFORT MASK

M-268020-21 (size S), M-268030-21 (size M), M-268040-21 (size L)

#### Packaging


# PACKAGING / LABEL BREEZE FACIAL COMFORT MASK

BREEZE NASAL+ MASK

M-268020-11 (size S), M-268030-11 (size M), M-268040-11 (size L)

## Definition of label stuck on the plastic bag

**BREEZE FACIAL Comfort Mask**  
Masque facial / Facial mask / Gesichtsmaske /  
Gezichtsmasker / Máscara facial / Maschera facciale

 0434			HSINER Co., Ltd. No.312, Jhongshan Rd., Shengang Dist. Taichung City 429
<b>REF</b> M-268020-21			MADE IN TAIWAN
<b>LOT</b> (10) 2013072513	<b>EC REP</b>		Mdi Europa GmbH Langenhagener Str. 71 30855 Langenhagen - GERMANY
 (11) 141017		<b>RX ONLY</b>	
 (17) 191017			

## PACKAGING / LABEL BREEZE NASAL COMFORT MASK

BREEZE NASAL COMFORT MASK

M-268020-11 (size S), M-268030-11 (size M), M-268040-11 (size L)

### Packaging



## PACKAGING / LABEL BREEZE NASAL COMFORT MASK

BREEZE NASAL+ MASK

M-268020-11 (size S), M-268030-11 (size M), M-268040-11 (size L)

### Definition of label stuck on the plastic bag

<b>BREEZE NASAL Comfort Mask S</b>			
Masque nasal / Nasal mask / Nasenmaske / Neusmasker / Máscara nasal / Maschera nasale			
 0434			HSINER Co., Ltd. No.312, Jhongshan Rd., Shengang Dist. Taichung City 429
 M-268020-12	MADE IN TAIWAN		
 (10) 2013072507		Mdi Europa GmbH Langenhagener Str. 71 30855 Langenhagen - GERMANY	
 (11) 141017			 -20°C / 60°C
 (17) 191017			

## PACKAGING / LABEL BREEZE NASAL MASK

BREEZE NASAL MASK

M-268020-10 (size S), M-268030-10 (size M), M-268040-10 (size L)

### Packaging



## PACKAGING / LABEL BREEZE NASAL MASK

BREEZE NASAL MASK

M-268020-10 (size S), M-268030-10 (size M), M-268040-10 (size L)

### Definition of label stuck on the plastic bag

<b>BREEZE NASAL Mask</b>			
Masque nasal / Nasal mask / Nasenmaske / Neusmasker / Máscara nasal / Maschera nasale			
 0434			HSINER Co., Ltd. No.312, Jhongshan Rd., Shengang Dist. Taichung City 429
 M-268020-10	MADE IN TAIWAN		
 (10) 2013072500		Mdi Europa GmbH Langenhagener Str. 71 30855 Langenhagen - GERMANY	
 (11) 141017		<b>RX ONLY</b>	 -20°C / 60°C
 (17) 191017			

## PACKAGING / LABEL BREEZE NASAL + MASK

BREEZE NASAL+ MASK

M-268020-11 (size S), M-268030-11 (size M), M-268040-11 (size L)

### Packaging


## PACKAGING / LABEL BREEZE NASAL + MASK

BREEZE NASAL+ MASK

M-268020-11 (size S), M-268030-11 (size M), M-268040-11 (size L)

### Definition of label stuck on the plastic bag

<b>BREEZE NASAL+ Mask</b>			
Masque nasal / Nasal mask / Nasenmaske / Neusmasker / Máscara nasal / Maschera nasale			
 0434			HSINER Co., Ltd. No.312, Jhongshan Rd., Shengang Dist. Taichung City 429
<b>REF</b> M-268020-11			MADE IN TAIWAN
<b>LOT</b> (10) 2013072503	<b>EC REP</b>		Mdi Europa GmbH Langenhagener Str. 71 30855 Langenhagen - GERMANY
 (11) 141017		<b>RX ONLY</b>	
 (17) 191017			

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
*ANNEX II excluding section 4 Directive 93/42/EEC concerning medical devices*

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*For class III devices, a EC design certificate is required*

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*Positive airways pressure devices for sleep disordered breathing treatment, humidifiers, software and associated accessories*

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*LNE/G-MED certifies that, on the basis of the results contained in the file referenced P116909, the quality system - for design, manufacturing, and final inspection - of medical devices listed here above complies with the requirements of the Directive 93/42/EEC, annex II excluding section 4*

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The validity of the certificate is subject to periodic or unexpected verification

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**Valable jusqu'au / Expiry date : March 15th, 2018 (included)**

**On behalf of the Deputy Director  
Cécile VAUGELADE  
G-MED Certification Division Manager**

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*Positive airways pressure devices for sleep disordered breathing treatment, humidifiers, software and associated accessories.*

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*The products covered by this certificate are listed on the list  
authenticated by LNE/G-MED on June 18<sup>th</sup>, 2015*

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*This certificate covers the following sites and activities*

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# EC DECLARATION OF CONFORMITY

For the following equipment:

Product Group: CPAP Mask and spare parts

Product Model: As Appendix

GMDN Code: 35174

The above product series are herewith confirmed to comply with the requirement set out in the Council Directive on the harmonization of the laws of member States concerning Medical Devices Directive (93/42/EEC-M5) described in Article 11.3 a) and Annex II.

This declaration is supported by the EC Certificate No. 47576-2009-CE-RGC-A 3.0 issued by notified body DNV (0434), and TCF No.TCF-03

For the evaluation regarding the products, the following standards were applied:

The following standards are applied for all our products:

EN 980:2008.

EN ISO 14971:2012.

EN ISO 10993-1:2009

EN ISO10993-5:2009.

EN ISO10993-10:2013.

EN 15986:2011

IEC 62366:2015

The following standards are applied to specific products as below:

Standards	Product
<u>BS EN ISO 17510-2:2009</u>	Sleep apnoea breathing therapy. Masks and application accessories
<u>BS EN ISO 5356-1:2015</u>	Anaesthetic and respiratory equipment. Conical connectors. Cones and sockets
<u>BS EN 20594-1:1994</u>	Conical fittings with a 6% (Luer) taper for syringes, needles and certain other medical equipment

Classification and rule of classification of MD: Class IIa (Annex IX, rule 2)

The following representative in Europe is responsible for this declaration:

mdiEuropa

Langenhagener Strasse71, 30855 Langenhagen, Germany

(Representative name & address)

Person responsible for making this declaration:

Manufacturer name; HSINER CO., LTD.

Manufacturer address: No.312, Jhongshan Rd., Shengang Dist., Taichung City 429, Taiwan

General Manager

AUG 08, 2016

(Position / title)

(Legal signature )

(Date)

## Appendix

HSINER REF	Description
M-268020-10	Breeze Nasal SIZE S
M-268030-10	Breeze Nasal SIZE M
M-268040-10	Breeze Nasal SIZE L
M-268020-11	Breeze Nasal+ SIZE S
M-268030-11	Breeze Nasal+ SIZE M
M-268040-11	Breeze Nasal+ SIZE L
M-268020-12	Breeze Nasal Comfort SIZE S
M-268030-12	Breeze Nasal Comfort SIZE M
M-268040-12	Breeze Nasal Comfort SIZE L
M-262040-00	Headgear Breeze Nasal, Nasal+&Facial+
M-262040-10	Headgear Breeze Comfort
M-429400-00	Silicone cushion Breeze Nasal S
M-429400-01	Silicone cushion Breeze Nasal M
M-429400-02	Silicone cushion Breeze Nasal L
M-429400-03	Silicone cushion Breeze Nasal+ S
M-429400-04	Silicone cushion Breeze Nasal+ M
M-429400-05	Silicone cushion Breeze Nasal+ L
M-429400-06	Silicone cushion Breeze Nasal Comfort S
M-429400-07	Silicone cushion Breeze Nasal Comfort M
M-429400-08	Silicone cushion Breeze Nasal Comfort L
M-267030-00	Elbow Breeze Nasal
M-267030-10	Elbow Breeze Nasal+&Comfort
M-367030-00	Elbow Breeze Nasal (x5)
M-367030-10	Elbow Breeze Nasal+ and Comfort (x5)
M-263030-10	Shell Breeze Nasal
M-263030-11	Shell Breeze Nasal+
M-263030-12	Shell Breeze Nasal Comfort
M-429410-00	Forehead Pad Breeze Nasal S
M-429410-01	Forehead Pad Breeze Nasal M
M-429410-02	Forehead Pad Breeze Nasal L
M-329410-00	Forehead Pad 3 sizes Breeze Nasal+
M-429410-10	Forehead Pad Nasal Comfort Mask
M-329411-00	Forehead Pad Breeze Nasal S (x3)
M-329411-01	Forehead Pad Breeze Nasal M (x3)
M-329411-02	Forehead Pad Breeze Nasal L (x3)
HSINER REF	Description

## Appendix

M-329411-03	Forehead Pad Breeze N+&F+ S (x3)
M-329411-04	Forehead Pad Breeze N+&F+ M (x3)
M-329411-05	Forehead Pad Breeze N+&F+ L (x3)
M-329411-10	Forehead Pad Breeze Nasal Comfort (x3)
M-268020-20	Breeze Facial+ SIZE S
M-268030-20	Breeze Facial+ SIZE M
M-268040-20	Breeze Facial+ SIZE L
M-268020-21	Breeze Facial Comfort SIZE S
M-268030-21	Breeze Facial Comfort SIZE M
M-268040-21	Breeze Facial Comfort SIZE L
M-429400-09	Silicone cushion Breeze Facial+&C S
M-429400-10	Silicone cushion Breeze Facial+&C M
M-429400-11	Silicone cushion Breeze Facial+&C L
M-263030-20	Shell Breeze Facial+ S
M-263030-21	Shell Breeze Facial+ M&L
M-263030-22	Shell Breeze Facial Comfort S
M-263030-24	Shell Breeze Facial Comfort M&L
M-268090-30	Breeze Pillows Mask
M-329400-02	Set of prongs for Breeze PM XS (x3)
M-329400-02	Set of prongs for Breeze PM S (x3)
M-329400-03	Set of prongs for Breeze PM M (x3)
M-329400-04	Set of prongs for Breeze PM L (x3)
M-262040-30	Headgear Breeze Pillows
M-429440-00	Clips Headgear Breeze N,N+&F+ (x10)
M-429440-01	Clips for Headgear Breeze Comfort (x10)
M-267030-30	Elbow Breeze wo leakage

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*Design, manufacturing, final inspection, commercialization and servicing of positive airways pressure devices for sleep disordered breathing treatment, humidifiers, software and associated accessories (see addendum).*

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**est conforme aux exigences des normes internationales**  
*complies with the requirements of the international standards*

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Renouvele le certificat 23443-3

**For the General Director**  
**Laurence DAGALLIER**  
**Deputy Director**

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*Design, manufacturing , final inspection, commercialization and servicing of positive airways pressure devices for sleep disordered breathing treatment, humidifiers, software and associated accessories.*

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## LETTER OF AUTHORIZATION & Distribution

We **SEFAM** authorize **MOBIAK S.A, 96-98 M. MPOTSARI Str, P.C: 73136 Chania, Crete, Greece,** to

1. Sell, Promote, and Distribute our Products.
2. Register our products to the National Organization for Medicines of Greece. :
  - DreamStar PAP Devices and accessories
  - Breeze masks products and accessories

The legal name and address of the manufacturer is:

SEFAM,  
144 avenue Charles de Gaulle,  
92200 Neuilly sur Seine  
FRANCE

The legal name and address of the manufacturing site:

SEFAM  
10 allée Pelletier Doisy  
54600 Villers-lès-Nancy  
France

We **SEFAM** declare that we sell and export our products, without any restriction, except Greece in the following three European Countries:

1. SPAIN
2. FRANCE
3. GERMANY

For and on behalf of the SEFAM company,

Pierrick HAAN

President