


FOLDING ANTI-SLIP WHEELCHAIR RAMP - 152CM

Κωδ: 0805403

- Length: 152 cm
- Inner Width of Each Part: 36cm
- Outer Width: 73cm
- Weight Tolerance: 270Kg



ROLLING RAMP

Κωδ: 0805401

- Length: 93cm
- Inner Width: 76cm
- Weight Tolerance: 270Kg



POWER WHEELCHAIR "VENERE"

Κωδ: 0808714

Venere is a Power Wheelchair with a Comfortable Seat, Made of Aluminium Alloy & with an Excellent Manufacturing Quality. It is Folding with Retractable Sides, Adjustable Footrest, and a Tilting Backrest for Easy Transport & Storage.



Technical Specifications

FRAME	Aluminium - Folding
COLOUR	Blue
SEAT	Fabric with reinforced blue-black cushion
ARMRESTS	Folding
FOOTREST	Detachable with Adjustable Length
REAR WHEELS	12" Solid
FRONT WHEELS	8" x 2" Solid (200x50mm)
SPEED	6-10Km/h
AUTONOMY	20-25Km
MAX TILT	10°
MOTOR	200W x 24 V x 2 Pieces / 1-Year Warranty
BATTERY	36Ah x 12V x 2 Pieces / 6-Month Warranty Ref: 0808529
CHARGER	24V x 4AMP / 1-Year Warranty
BRAKES	Magnetic

- Right- or Left-Handed Controller
- Anti-tilt Casters
- Safety Belt
- Tilting Backrest

0808714	46cm	45cm	41cm	45cm	62cm	106cm	93cm	51cm	12" 8"	110kg	49kg

"HERMES II" POWER WHEELCHAIR

Κωδ: 0811315

Hermes is a Power Wheelchair with a Comfortable Seat, Folding, with Adjustable Footrest and a Tilting Backrest for Easy Transport & Storage.



- Right- or Left-Handed Controller
- Anti-tilt Casters
- Tilting Backrest

Technical Specifications

FRAME	Steel - Folding
COLOUR	Dark Grey
SEAT	Fabric with Reinforced Cushion
ARMRESTS	Folding
FOOTREST	Detachable with Adjustable Length
REAR WHEELS	12.5" Solid
FRONT WHEELS	8" x 2" Solid (200x50mm)
SPEED	7Km/h
AUTONOMY	15Km
MAX TILT	10°
MOTOR	250W x 24 V x 2 Pieces / 1-Year Warranty
BATTERY	20Ah x 12V x 2 Pieces / 6-Month Warranty Ref: 0804924
CHARGER	24V x 3AMP / 1-Year Warranty
BRAKES	Magnetic

0811315	120kg	45cm	44cm	44cm	63cm	108cm	97cm	53cm	12,5" 8"	120kg	35kg	

"TITAN" POWER WHEELCHAIR

Κωδ: 0811317

Titan is an Exceptionally Comfortable Power Wheelchair with Many Different Adjustments that allow Absolute Independence of Movement!



Technical Specifications

FRAME	Aluminium - Folding
COLOUR	Silver
SEAT	Removable Cushion
ARMRESTS	Detachable - Height Adjustable
FOOTREST	Detachable with Adjustable Length
REAR WHEELS	13.5" Solid
FRONT WHEELS	9" x 2" Solid (230x50mm)
SPEED	8-10Km/h
AUTONOMY	30-36Km
MAX TILT	10°
MOTOR	320W x 24 V x 2 Pieces / 1-Year Warranty
BATTERY	40Ah x 12V x 2 Pieces / 6-Month Warranty Ref: 0804925
CHARGER	24V x 5AMP / 1-Year Warranty
BRAKES	Magnetic

- Right- or Left-Handed Controller
- Anti-tilt Casters
- Safety Belt
- Tilting Backrest

0811317	46cm	44cm	50cm	40cm	65cm	110cm	100cm	55cm	13.5" 9"	130kg	67kg

"VOYAGER" POWER WHEELCHAIR` FOLDABLE

Κωδ: 0811316

Voyager is a Power Wheelchair, Small in Size, Made of Aluminium Alloy, Folding, Ideal to Go Everywhere & Easy to Store.



Technical Specifications

FRAME	Aluminium - Folding
COLOUR	Silver
SEAT	Removable Cushion
ARMRESTS	Folding
FOOTREST	Folding
REAR WHEELS	12" Solid
FRONT WHEELS	8" x 2" Solid (200x50mm)
SPEED	7Km/h
AUTONOMY	15-20Km
MAX TILT	10°
MOTOR	200W x 24 V x 2 Pieces / 1-Year Warranty
BATTERY	20Ah x 24V x 1 Piece / 6-Month Warranty Ref: 0804923
CHARGER	24V x 3AMP / 1-Year Warranty
BRAKES	Magnetic

- Right- or Left-Handed Controller
- Anti-tilt Casters
- Safety Belt
- Lithium Batteries

0811316	48cm	43cm	43cm	45cm	65cm	110cm	105cm	55cm	12 8"	120kg	38kg