



Code: 0806270

"MINOA" Anti-Decubitus Bubble Mattress with Electric Air Pump

- Adjustable Electric Pressure Control according with the user's weight requirements
- Maximum Patient Weight: 130Kg
- Mattress Dimensions (LxWxH): 200 x 90 x 7cm
- Mattress can be folded under the foam mattress with 2 flaps
- Pump Dimensions (LxWxH): 23.5 x 11.7 x 8cm
- Mattress Warranty: 6 Months
- Pump Warranty: 2 Years
- Product Weight: 3Kg



Στρώμα Κατακλίσεων
με Αεραντλία "MINOA"
Anti - Decubitus Air Mattress
With Pump "MINOA"
Dušek "MINOA"

Matelas à air de
prévention-escarre "MINOA"
Anti-decubitus
luchtbed "MINOA"
Die Luftmatratze
"MINOA"

REF: 0806270
ΚΩΔ: 0806270

MOBIAM S.p.A.
KAROLINA APOSTOLIOU, KAROLINA APOSTOLIOU, KAROLINA APOSTOLIOU, KAROLINA APOSTOLIOU
T: +30 26110 60221 | WEB: www.mobiak.com

MD CE MOBIAM



Code: 0806276

"PREMIUM OCEAN" Anti-Decubitus Bubble Mattress with Electric Air Pump

- Adjustable Electric Pressure Control, according with the User's Weight: 30 - 130kg
- Static & Alternating modes
- Maximum Patient Weight: 130Kg
- Mattress Dimensions (LxWxH): 200 x 90 x 7cm
- Mattress can be folded under the foam mattress with 2 flaps
- Pump Dimensions (LxW): 26 x 12,5 x 8,5cm
- Mattress Warranty: 6 Months
- Pump Warranty: 2 Years
- Product Weight: 3,3Kg

