



Code: 0806380

Gem Nebulizer

- The Device has been Designed for Easy use, Transport and it's Recommended for Nebulization of Antibiotics and Bronchodilators
- Fast Nebulization, Saving up to 50% Time thanks to the High-flow Nebulization Chamber (Ref: 0148006)
- Includes a High-flow Mask for Adults & Children, an Accessory for Nose Use, a Highflow Nebulization Chamber and a Spare Filter



Technical Specifications

Max Pressure	2.5 bar
Operating pressure	1.1 bar
Max Flow	14Lt /min
Operating Air Flow	4.8Lt / min
Nebulization Flow with 4 ml Solution	0.35 ml/min
Mean Aerodynamic Mass Diameter	2.44 μm (MMAD)
Inhaled Ratio (<5μm)	> 80%
Supply	230 V, 50Hz
Power Consumption	170 VA
Operating Conditions	Intermittent (20 min ON - 40 min OFF)
Noise Level	60 dB
Device Weight	1.4Kg
Device Dimensions	13 x 8.9 x 14.6cm
Warranty	5 years

Package Contents

