

Auto-CPAP iBreeze Series 20A

Ref: 0803495



Your Auto CPAP Package include the following items:

- Main Device
- Power Adapter and Cables
- Flexible Air Tube
- User Manual
- Travel Bag
- Filter (2pcs)
- SD Card

iBreeze Series is the brand new generation of Sleep & Respiration Therapy Platforms developed by Resvent. It is more than a simple device, a systematic solution of Respiratory Disorders which includes Diagnosis, Titration and Treatment.

Ergonomic design, intuitive user interface and Advanced Ventilation Algorithm along with ResAssist (a Resvent's Cloud based patient data Management System), will bring users the most brilliant experience.

Technical Specifications

Modes Of Operation	CPAP, Auto-CPAP
Therapy Pressure	4-20cm H ₂ O
Ramp - Time	0-60min, 5min increment
Autoramp	Yes
Humidifier	Auto or 0-8 levels of adjustment
Sound pressure level	36dBA
Screen Size	3.5"
Dimensions (length*width*height)	238 x 178 x 128 mm
Weight	1,6kg
Warranty	2 Years



BiLevel

Ref: 0803496

Series 25A

Ref: 0803497

Series 25ST

Ref: 0803498

Series 25STA

MOBIAK®
MEDICAL • FIREFIGHTING • GASES est. 1977



Integrated Humidifier

Intelligent compact design with integrated humidifier brings more convenience in your daily treatment at home and maximizes your confidence during your travels.

Ergonomics

Innovative inclined front-panel's design, makes it easy for operation while lying in bed as well as sitting up on bedside.

Touch screen

With intuitive 5.0" touch screen display and knob navigation, there are no barriers between seeing and using.

Intelligence Pressure Release

Unique IPR (Intelligence Pressure Release) algorithm during inspiration & expiration offers comfortable experience to patients during ventilation therapy. Three IPR levels are available to meet your personal requirements.

REG

Innovative Algorithm by Resvent can continuously determine and monitor respiratory events (AHI, OAI, CAI, HI, SNI, FL, RERA, PB, etc), analyze the physical status of patients and guarantee the delivery of the optimal therapy pressure every time.

Auto-humidifier

During respiration treatment, water condensation in the breathing circuit affects the therapy and even threatens the safety of the patients. In order to avoid this, Resvent's Unique Intelligent humidifier automatically adjusts the degree of the humidifying according to the environment's temperature & humidity as to achieve zero condensation in the breathing circuit and/or to the mask.



BiLevel

Ref: 0803496

Series 25A

Ref: 0803497

Series 25ST

Ref: 0803498

Series 25STA

MOBIAK®
MEDICAL • FIREFIGHTING • GASES est. 1977



Your BiLevel Package include the following items:

- Main Device
- Power Adapter and Cables
- Flexible Air Tube
- User Manual
- Travel Bag
- Filter (2pcs)
- SD Card

Technical Specifications

Item	Value		
	25A	25ST	25STA
Modes Of Operation	CPAP ,S ,Auto S	CPAP ,S ,T, S/T	CPAP ,S , Auto S,T, S/T, Auto S/T
Therapy Pressure	4-25cm H ₂ O		
Ramp- Time	0-60min, 5min increment		
Humidifier	Auto or 0-8 levels of adjustment		
Sound pressure level	36dBA		
Screen	5", Touch		
Dimensions (LxWxH)	238 x 178 x 128 mm		
Weight	1,65 kg		
Warranty	2 Years		

All iBreeze Series devices offer comprehensive and adaptive connectivity options for reviewing the therapy data. Patients and Doctors can easily be connected to Resvent's Cloud, Data Server using both Wi-Fi and/or GSM Module. There is also a traditional SD card available for storage of the therapy data.

Resvent AUTO-CPAP iBreeze 20A

Ref: 0803494



Your Auto CPAP Package include the following items:

- Main Device
- Power Adapter and Cables
- Flexible Air Tube
- User Manual
- Travel Bag
- Filter
- SD Card

The iBreeze 20 A is one of the smallest portable Auto Positive Airway Pressure (Auto CPAP) systems in the world. With its help you can overcome all the obstacles between healing and freedom. Regardless of its small size, it can withstand even the most demanding treatments.

Technical Specifications

Modes Of Operation	CPAP, Auto-CPAP
Therapy Pressure	4-20cm H ₂ O
Ramp- Time	0-60min, 5min increment
Power Supply	100-240VAC, 50/60Hz
Battery Duration	(3200mAh) - 6h // (6600mAh) - 12h
Operation Temperature	5 ~35 °C
Sound Level	36 dB(A)
Dimensions (LxWxH)	151 x 151 x 72 cm
Weight (Kgr)	0.68 Kgr
Warranty	2 Years

