



Oxygen Concentrator 5lt Irene

Code: 0808194

Technical Specifications

Oxygen Purity (%)	93± 3
Oxygen Flow (LPM)	1-5
Max Output Pressure (bar)	0,5
Input Voltage (VAC)	220V± 10%, 50Hz± 1
Power Consumption (W)	480
Sound Level dB(A)	≤40
Dimensions (cm)	L:39 x W:33,7 x H:62
Weight (Kg)	24

- Large Compressor for More Strength and Durability

Package Contents

- Nasal Cannula
- Humidifier
- Humidifier Adaptor & Tube
- Cabin Filter
- Inlet Filter
- User Manual
- Power Supply Cable



Oxygen Concentrator 10lt FORCE

Code: 0810360

Technical Specifications

Model	FORCE
Power (W)	635
Voltage/Frequency	220V± 22V 50HZ ± 1 HZ
Flow rate (Lt/min)	0-10
Purity at 5 l/min	93+3%
Purity at 8 l/min	90+3%
Purity at 10 l/min	87+3%
Output Pressure (Mpa)	0.05±0.005
Noise Level dB(A)	≤55
Classification & Type	Class II Type A
Unit Weight (Kg)	27
Dimensions (mm)	381x347x689

Package Contents

- Nasal Cannula
- Humidifier
- Cabin Filter
- Inlet Filter
- User Manual
- Power Supply Cable



Oxygen Concentrator 5lt 8F-5A

Code: 0803300

Technical Specifications

Oxygen Purity (%)	93± 3
Oxygen Flow (LPM)	0.5-5
Max Output Pressure (bar)	0,7
Input Voltage (VAC)	220V± 10%, 50Hz± 1
Power Consumption (W)	400
Sound Level dB(A)	53
Dimensions (cm)	L:39 x W:24,5 x H:50
Weight (Kg)	16,6

Package Contents

- Remote Control
- Power Cord (EU Type)
- Air Intake Filter
- User Manual



Oxygen Concentrator 5lt AR-5

Code: 0804199

Technical Specifications

Oxygen Purity (%)	93± 3
Oxygen Flow (LPM)	1-5
Max Output Pressure (bar)	0,5
Input Voltage (VAC)	220V± 10%, 50Hz± 1
Power Consumption (W)	390
Sound Level dB(A)	≤45
Dimensions (cm)	L:31,4 x W:24 x H:51,5
Weight (Kg)	13,9

Package Contents

- Nasal Cannula
- Humidifier
- Humidifier Adapter and Tube
- Cabin Filter
- Air Intake Filter
- Power Cord (EU Type)
- User Manual

