

FreeStyle® Comfort®

The perfect fit



Key features and benefits



Optimal Comfort

Let your patients experience the comfort of the body-skimming, curved silhouette design. The FreeStyle Comfort portable oxygen concentrator is lightweight and compact.




Smart O2 delivery

Unique clinical features including UltraSense®, autoDOSE, and myCAIRE™ help to ensure your patients will get the oxygen they need, when they need it.



Proper O2 Saturation

With 5 settings and up to 1050 mL of oxygen per minute, the FreeStyle Comfort serves a wide range of patient needs. The unit delivers up to 95% pure oxygen to further ensure proper saturation.



Proven Quality

Robust design features, including the sieve bed technology and critical components, allow for minimum maintenance and long life. CAIRE offers one of the best warranties available, with a 3-year warranty on the device and sieve beds.

Specifications



FLOW RATE

Pulse Dose 1-5
O2 Output Settings:
Setting 1-210 mL/min
Setting 2-420 mL/min
Setting 3-630 mL/min
Setting 4-840 mL/min
Setting 5-1050 mL/min

DIMENSIONS

10.0 x 7.3 x 3.1 in
(25.4 x 18.5 x 7.9 cm)
8-cell battery pack
11.0 x 7.3 x 3.1 in
(27.9 x 18.5 x 7.9 cm)
16-cell battery pack

WEIGHT

5 lbs (2.3 kg)
8-cell battery pack
6 lbs (2.7 kg)
16-cell battery pack

SOUND LEVEL

39.9 dB(A) @ Setting 2

BATTERY DURATION

8-cell battery
up to 8 hours
16-cell battery
up to 16 hours

POWER

AC Power: 100-240
VAC, 50-60Hz
DC Power: 11-18 VDC
(10 max amp)
Battery: Lithium Ion

WARRANTY

3 years
(Includes sieve beds)
1 year on battery and
power supplies