



COMPRESSION THERAPY



The Air Compression Therapy System provides a High-Quality Lymphatic Decongestion. It also Helps Blood and Lymph Circulation, Preventing Deep Vein Thrombosis (DVT) and Provides Swelling Relief. It Crucial Contributes to Blood and Lymph Circulation in Arms and Legs by Enhancing Muscle Tissue Pressure to Remove Swelling or Pain. It is Used for Prevention and Treatment of Various Conditions Related to Arterial Blood Circulation. It is Used by persons who suffer from Diabetes or Post-Thrombotic Syndrome.

Prescription-Only – Administered by Qualified Technicians or Nursing Staff.

One or two Pads can be Connected to the same Compressor (if it is suggested by the Doctor).



Leg Compression Pad (4 Chambers)

Ref: 0806213

Ref: 0806214

Ref: 0806225



Compressor, Digital (4 Ports)

Ref: 0806211


- Adjustable Pressure (Pressure range 20 - 250 mmHg)
- User-Selected Compression Therapy Duration
- 6 Programmes Available



Compressor, Analogue (4 Ports)

Ref: 0806210


- Adjustable Pressure (Pressure range: 30 - 240mmHg)
- Selection of Compression Therapy Duration
- Purchase of Suitable Pad



Arm Compression Pad (4 Chambers)

Ref: 0806216

Ref: 0806217



2-Pad Connector (4 Chambers)

Ref: 0806221



Compression System Carrying Bag

Ref: 0806220



Back Compression Pad (4 Chambers)

Ref: 0806213

Ref: 0806214



www.mobiakcare.com

Value in safety, Invest in Quality!