



CYAN COMMUNE WHEELCHAIR

Ref: 0805453

SPECIFICATIONS

High-quality wheeled commode chair combining easy-to-clean upholstery with a reliable and sturdy frame.

- Frame: Steel
- Adjustable - Folding Armrests
- Detachable, Adjustable Length & Rotating Footrests



One-handed bucket operation



PVC Seat



4" Solid Wheels with brakes for safety



0805453	44,5cm	43cm	43cm	46cm	55cm	77cm	97cm	53cm	4" 4"	136Kg	11Kg	

www.mobiakcare.com



Value in safety, invest in Quality!