

Reinforced Foldable walkers with ergonomic handles, suitable for all users with weight until **136kg**.

New Product



Foldable Walker 136KG

Ref: 0809786

SPECIFICATIONS

- Frame: Foldable with one button
- Color: Silver – Black
- Light-weight: 2,4 kg
- Height-Adjustable to cover every adult user, height demand
- Anti-slip Ferrule, Ref: 0809788



58,5cm	42,5cm	47cm	75-93cm	136Kg

Foldable Walker with “Step-By-Step” Motion 136KG

Ref. 0809787

SPECIFICATIONS

- Frame: Foldable with one button
- With Reciprocating Motion, with stabilization possibility
- Color: Silver – Black
- Light-weight: 2,8 kg
- Height-Adjustable to cover every adult user, height demand
- Anti-slip Ferrule, Ref: 0809788

New Product



52cm	44,5cm	47cm	75-93cm	136Kg