

Aluminium Wheelchair AZURA

Code: 0811982

Code: 0811983



Specifications

- Aluminum Frame
- Frame Color Blue
- Reinforced Seat from Fabric
- Detachable Armrests with Adjustable Length
- Swing Away, Reclining & Adjustable Height Footrests
- Rear Wheels 16" Solid
- Front Wheels 6" x 1" Solid (152x25mm)
- Max Load 100kg
- 45° Reclining Seat
- 90° to 165° Backrest Tilt
- Safety Cross Belt
- Height-adjustable, Anti-tilt Support Casters

A Wheelchair suitable for people with paraplegia or quadriplegia (cerebral palsy), provides comfort for various postures, this is favored by the possibility of adjusting the headrest, side panels and stern protectors. The headrest, the side panels, the sternum supports, and the foot divider are lined with soft material for greater comfort and the cross-belt offers bigger safety. It is foldable and facilitate the user for his/her transportation by car.

| 0811982 | 43-49cm | 42cm | 40cm | 45cm* | 62-64cm | 102-114cm** | 104-117cm*** | 50cm | 16" 6" | 100kg | 21,5kg |
|----------------|---------|------|------|-------|---------|-------------|--------------|------|--------|-------|--------|
| 0811983 | 38-44cm | 36cm | 40cm | 45cm* | 58-60cm | 102-114cm** | 104-117cm*** | 50cm | 16" 6" | 100kg | 21kg |

*without the headrest

**with the footrests in the lower and fully reclining position

***with the headrest in the lower and higher position